

# Using EquiChek™-SAA for a diagnostic follow up exam



EquiChek™-SAA can be conveniently run while you perform a follow up exam. No sample processing is needed. With just a pin prick of blood you can detect and monitor acute inflammatory conditions in less than 10 minutes.

## EquiChek™-SAA + Physical Exam



4 lines means no inflammation  
(Recovery)



3 Lines means mild inflammation



2 Lines means moderate inflammation



1 Line means severe inflammation, possible infection\*

Recovery from earlier illness

Recovering from earlier illness

Recovery started monitor for 3 - 7 days

Ongoing active inflammatory condition

Confirmation of recovery

Continue to monitor

If repeatedly 2 lines suspect internal disease or infection

Consider other symptoms and recheck in two days

If no change in levels reconsider diagnosis and treatment options

\*A definitive clinical diagnosis should only be made after all clinical and laboratory findings have been considered. Information is provided for user understanding and does not constitute veterinary advice.

Contact: Equine Partners America, LLC  
Tel: +1(800) 752-8538, Cell: +1(770) 530-5670  
Email: moira@equinepartnersamerica.com  
Web: www.EquinePartnersAmerica.com

