

How to use OmniChek™-SAA in a diagnostic workup



OmniChek™-SAA is a simple to run rapid test for inflammation and infection that can conveniently be performed while you conduct a physical exam. No sample processing is needed. With just a pin prick of blood you can detect and monitor acute inflammatory conditions in less than 10 minutes.

OmniChek™-SAA



3 lines

Non
Inflammatory/
Non Infectious

Recovery
from earlier
illness

Proceed with
examination
having ruled
out these
conditions

Confirmation
of recovery to
discontinue
medication



2 lines

Mild/Moderate
inflammation in most
species

If a first
measure,
consider case
history

If repeatedly 2
lines suspect
internal
disease or
neoplasia



1 line

Severe inflammation,
possible infection*

Physical exam,
history and
further testing
to confirm
cause

Consider
commencing
treatment
based on
evidence

*A definitive clinical diagnosis should only be made after all clinical and laboratory findings have been considered. Information is provided for user understanding and does not constitute veterinary advice.

Contact: Equine Partners America, LLC
Tel: +1 (800) 752-8538, **Cell:** +1 (770) 530-5670
Email: moira@equinepartnersamerica.com
Web: www.EquinePartnersAmerica.com

