

# Using OmniChek™-SAA for a diagnostic follow up exam



OmniChek™-SAA can be conveniently run while you perform a follow up exam. No sample processing is needed. With just a pin prick of blood you can detect and monitor acute inflammatory conditions in less than 10 minutes.

## OmniChek™-SAA + Physical Exam



3 lines means no inflammation  
(Recovery)



2 lines means mild/moderate inflammation



1 line means severe inflammation, possible infection\*

Recovery from earlier illness

Recovery started, monitor for in 3 - 7 days

If this is a persistent result over more than 5 days consider treatment options. Ongoing active inflammatory condition

Confirmation of recovery

If repeatedly 2 lines suspect internal disease or neoplasia

Consider other symptoms and recheck in 2 days

If no change in levels, reconsider diagnosis and treatment options

\*A definitive clinical diagnosis should only be made after all clinical and laboratory findings have been considered. Information is provided for user understanding and does not constitute veterinary advice.

**Contact:** Equine Partners America, LLC  
**Tel:** +1 (800) 752-8538, **Cell:** +1 (770) 530-5670  
**Email:** moira@equinepartnersamerica.com  
**Web:** www.EquinePartnersAmerica.com

